---

---

# SKED<sup>®</sup>

## THE COMPLETE RESCUE SYSTEM

MANUAL

# SKED® SAVES LIVES

Congratulations on your purchase of the Sked® Rescue System. You will find it a valuable tool that will enable rescues to be performed safer, faster, and easier. The Sked® Stretcher is an entirely different concept in stretcher design. Therefore, potential rescuers should practice and familiarize themselves with the Sked® prior to using it on an actual rescue.

The Sked® Stretcher provides excellent patient support and protection. However, it is NOT designed as a spinal immobilization device. If a spinal injury is suspected, secure the patient to an approved spinal immobilization device prior to placing him/her in the Sked®. The Sked® Stretcher will accommodate long and short backboards, scoop stretchers, Oregon Spine Splint II™ and most other immobilization equipment. A backboard must be used in conjunction with the Sked® Stretcher on patients who have sustained injuries to the shoulder area.

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# TO UNROLL SKED® STRETCHER AND LAY FLAT

\*STEEL BUCKLES


\*COBRA BUCKLES

2. UNFASTEN RETAINER STRAP. STEP ON FOOT END OF SKED® AND UNROLL COMPLETELY TO OPPOSITE END.

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1. REMOVE SKED® FROM PACK AND PLACE ON GROUND.

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3. BEND THE SKED® OVER YOUR FOREARM AND BACK ROLL. REPEAT WITH OPPOSITE END OF STRETCHER. IN COLD WEATHER, IT MAY BE NECESSARY TO REVERSE ROLL IT TWICE.

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# PLACE A PATIENT ON THE SKED®: PULL METHOD (FEET FIRST)

1. PRIOR TO DRAGGING A PATIENT ONTO THE SKED® "FEET FIRST". CURL THE HEAD END DOWNWARD TO FORM A RAMP AND TO PREVENT SNAGGING PATIENT'S CLOTHING.

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2. TUCK THE DRAG WEBBING UNDER THE SKED®. GRASP PATIENT'S LEGS AND DRAG PATIENT WHILE KEEPING LEGS LOW.

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3. CONTINUE DRAGGING UNTIL HIS UPPER CHEST IS EVEN WITH THE UPPER CROSS STRAP.

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4. SECURE ALL BUCKLES AND ADJUST ALL STRAPS IN PREPARATION FOR TRANSPORT.

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# PLACE A PATIENT ON THE SKED®: PULL METHOD (HEAD FIRST)

1. CURL FOOT END OF SKED® DOWNWARD TO FORM A RAMP AND PULL FOOT END STRAPS TO THE SIDE.

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
# FASTEN THE STRAPS: COBRA QUICK RELEASE BUCKLES

1. PRIOR TO LOADING THE SKED® RELEASE ALL COBRA BUCKLES AND PULL STRAPS TO ONE SIDE.

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2. LIFT SIDES OF SKED® AND FASTEN THE FOUR CROSS STRAPS TO BUCKLES DIRECTLY OPPOSITE TO THE STRAPS AND PULL THEM TIGHT.

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3. FASTEN FOOT END COBRA BUCKLES AND ADJUST STRAPS.

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
4. DO NOT CRISS-CROSS STRAPS ON THE SKED®. IT WILL CAUSE THE SKED® TO BEND IN THE MIDDLE.

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# FASTEN THE STRAPS: STEEL BUCKLES

1. PRIOR TO LOADING THE SKED® RELEASE ALL STEEL BUCKLES AND PULL STRAPS TO ONE SIDE. TO RELEASE, LIFT THE LEADING END OF THE BUCKLE TO LOOSEN WEBBING. PULL THE WEBBING OUT OF THE BUCKLE.

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2. FASTEN THE 4 CROSS STRAPS TO BUCKLES DIRECTLY OPPOSITE TO THE STRAPS. TO CONNECT, SEE STEP #3 TO LACE BUCKLES, THEN TIGHTEN STRAPS.

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4. FEED THE STRAPS THROUGH UNUSED GROMMETS AT FOOT END OF SKED AND FASTEN TO BUCKLES THEN ADJUST STRAPS.

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# HORIZONTAL LIFT/DESCENT

1. INSERT ONE END OF HEAD STRAP THROUGH ANGLED SLOT AT HEAD END OF SKED®.

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2. BRING STRAP UNDER SKED® AND THROUGH ANGLED SLOT ON OPPOSITE SIDE OF SKED®.

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3. EQUALIZE STRAP. REPEAT PROCEDURE WITH OTHER STRAP AT FOOT END OF SKED®.

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
4. EQUALIZE ALL FOUR STRAPS AND SECURE ENDS WITH LARGE STEEL LOCKING CARABINER.

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# VERTICAL LIFT/DESCENT

1. PRIOR TO HOISTING. THE HEAD END MUST BE CURLED OVER THE PATIENT'S HEAD AND THE DRAG HANDLE WEBBING TIED TO THE SECOND CROSS STRAP TO PROTECT PATIENT'S HEAD FROM FALLING OBJECTS.

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3. FOR EXTREMELY SMALL OPENINGS, LEAVE THE HEAD END EXTENDED. START THREADING THE ROPE THROUGH THE GROMMETS BY THE PATIENT'S SHOULDERS.

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
2. INSERT VERTICAL LIFT SLING ROPE THROUGH THE GROMMETS ABOVE PATIENT'S HEAD.

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4. CONTINUE THREADING THE ROPE ON EACH SIDE OF THE SKED® AS SHOWN IN PICTURE ABOVE.

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# VERTICAL LIFT/DESCENT

5. THREAD THE ROPE THROUGH THE BOTTOM GROMMETS FROM THE INSIDE-OUT AS SHOWN ABOVE.

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6. TIE THE ROPE IN A SQUARE KNOT AT THE LOWER END OF THE SKED®.

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7. PASS THE ENDS OF THE ROPE THROUGH THE LOWER CARRY HANDLES FROM THE OUTSIDE-IN AND TIE ANOTHER SQUARE KNOT.

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8. TIE AN OVERHAND KNOT ON EACH SIDE OF THE SQUARE KNOT FOR SAFETY.

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# BACKPACK/TOWING HARNESS

 1. ATTACH ONE END OF THE TOW STRAP TO THE BACKPACK. BACKPACK IS USED AS TOWING HARNESS.

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2. ATTACH THE OTHER END OF THE TOW STRAP TO SKED® DRAG WEBBING. THIS ALLOWS "HANDS FREE" DRAGGING OF THE SKED®.

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## ACCESSORIES: TOW STRAPS/REMOVABLE WEBBING HANDLES

 3. TOW STRAP IS ATTACHED TO CARRY HANDLES ENABLING 3 OR 4 RESCUERS TO DRAG THE SKED®.

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4. SET OF 4 REMOVABLE WEBBING HANDLES ENABLES UP TO 8 RESCUERS TO CARRY SKED®. INSERT THRU UNUSED GROMMETS ON SIDES OF SKED®.

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# ROLLING UP THE SKED® STRETCHER

1. LAY STRETCHER OUT FLAT WITH ALL CROSS STRAPS CONNECTED. \*\*STEEL DOUBLE BAR BUCKLES: PLACE RETAINER STRAP WITH INSIDE OF BUCKLE FACING UP UNDER FOOT END OF SKED®. (SEE PHOTO IN STEP 4)

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2. STARTING AT THE HEAD END, ROLL THE SKED® UP AS TIGHT AS POSSIBLE. CONTINUE TO ROLL THE SKED® USING KNEE TO KEEP THE ROLLED STRETCHER SMALL.

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3. WHILE HOLDING YOUR KNEE ON THE ROLLED SKED®, FASTEN THE COBRA BUCKLES OF THE FOOT END STRAPS TO HOLD THE SKED® IN A TIGHT ROLL.

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4. WHEN USING THE SKED® WITH STEEL DOUBLE BAR BUCKLES, FASTEN THE PRE-PLACED RETAINER STRAP TO BUCKLE AND PLACE SKED® STRETCHER IN BACKPACK.

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# HORIZONTAL LIFT WITH VERTICAL LIFT SLING

1. HOLD THE PRE-TIED FIGURE 8 KNOT IN THE CENTER OF THE ROPE AND PASS THE ENDS OF THE ROPE THROUGH THE ANGLED LIFT SLING SLOTS FROM THE INSIDE-OUT.

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3. TIE FIGURE 8 KNOTS IN BOTH ENDS OF THE ROPE LEAVING ABOUT 2.5 FEET EXTENDING UP FROM THE FOOT END SLOTS AND SHORTER AT THE HEAD END. (SEE PHOTO IN STEP 4)

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2. PASS THE ENDS OF THE ROPE UNDER THE SKED® AND THROUGH THE ANGLED LIFT SLING SLOTS ON THE OPPOSITE SIDE FROM THE OUTSIDE-IN.

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4. THE ROPE AT THE HEAD END MUST BE A BIT SHORTER THAN THE FOOT END TO PROVIDE A SLIGHTLY "HEAD UP" ATTITUDE FOR PATIENT'S SAFETY AND COMFORT.

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## What does “SKED” stand for?

Though we spell it with capital letters, SKED® is not an acronym. It came from fusing two words: “Sled” and “Skid.” The early idea behind the product was that it was a SLED that SKIDDED across all types of terrain. Initially, the SKED® was used as a game carrier to tote wild game back to a hunter’s camp. It has since morphed into the life-saving device that’s used today, worldwide. Like many iconic products that are the first of its kind in the market, the brand has come to epitomize the product itself.

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## When does the SKED stretcher need to be replaced?

We have a max time limit off 10 years, when the SKED® needs to be replaced because different rescue professionals use it in different ways for different lengths of time. If there is any damage to your stretcher, please reevaluate if you should be using it. If you have questions do not hesitate to call us.

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## What should I do if my SKED® is missing a grommet?

SKEDCO uses two grommet sizes that can be found at most local hardware stores. There are two different sizes: the smaller grommet is a size #3, and the larger is a size #5.

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## What should I do if the straps on my SKED® are starting to fray?

We have replaceable straps – the SK-208C. The SKED-EVAC® Aluminum Side Release Buckle Conversion Kit features Austrian-made Skedco/Cobra side release buckles. These buckles are dependable and strong enough to meet the needs of securing a patient into a SKED® stretcher without the danger of breaking the buckle or accidental release.

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# SKED® FACT SHEET

THE SK-200-OR SKED® BASIC RESCUE SYSTEM – INTERNATIONAL ORANGE COMES WITH:

Available with Flotation and the Oregon Spine Splint II (OSS) options

Although the SKED system works very well for virtually all rescue situations, it excels in confined space, urban, military and industrial rescue.

# SKED® FACT SHEET

## SIZE

Rolled in cordura pack – 23 cm diameter x 91cm length

Laid out flat – 91cm x 244cm

WEIGHT 9kg. complete with accessories listed below.

MATERIAL Low density E-Z glide polyethylene plastic.

\*\* Nylon webbing

\*\* Solid brass grommets

\*\* Steel buckles

TEMPERATURE RANGE Begins to melt 230 degrees C to 260 degrees C

\*\* Usability to - 84 degrees C without becoming brittle.

HORIZONTAL LIFT SLINGS Webbing tensile strength 45 KN.

VERTICAL LIFT SLING Rope tensile strength 31KN.

SYSTEM COMPONENTS SKED STRETCHER - BASIC RESCUE SYSTEM INCLUDES:

\*\* SKED Stretcher

\*\* Cordura Back Pack/Towing Harness

\*\* Horizontal Lift Slings

\*\* Vertical Lift Sling (PMI E-Z Bend 9mm rope)

\*\* Large steel locking carabiner (SMC)

\*\* Tow Strap

\*\* Set of 4 removable webbing handles

Available with Flotation and the Oregon Spine Splint II (OSS) options

# SKED® MAINTENANCE

## INSPECTION FORM VISUAL CHECKING OF THE SKED STRETCHER

TEXTILE PART : Condition of webbing (cuts, wear, burns) Condition of stitching (broken, pulled, or worn threads)  
 METALLIC PART : Condition of anchor-rings (deformation, marks, cracks, wear, corrosion) Condition of the adjustment buckles deformation, marks, cracks, wear, corrosion)  
 PLASTIC PART: : Plastic sheet (deformation, marks, cracks, wear, corrosion)

### OPERATIONAL CHECK G TM TR R C

G: Good / TM: To Monitor / TR: To Repair / R: Reject / C: Comment

Date of inspection :.....Date of next inspection :.....

Inspected by.....Company name:.....

Historical Check:.....

Year of manufacturing:.....

Date of Purchase:.....

Ind. serial nr:.....

| PRODUCT           | G | TM | TR | R | C |
|-------------------|---|----|----|---|---|
| BODY STRETCHER    |   |    |    |   |   |
| HORIZONTAL SLINGS |   |    |    |   |   |
| VERTICAL ROPE     |   |    |    |   |   |
| KARABINER         |   |    |    |   |   |
| EXTRA HANDLES     |   |    |    |   |   |
| TOW STRAP         |   |    |    |   |   |
| GROMMETS          |   |    |    |   |   |

COMMENTS.....

VERDICT: The product is fit ..... to remain in service. The product is unfit ..... to remain in service  
 In addition to routine checks for each use, Stretchers should regularly undergo a detailed inspection by a competent person. Skedco Europe recommends an inspection every 12 months and after any exceptional event in the life of the product. Inspections that do not come into any of the categories listed above, any of which would require the systematic rejection or the component, namely:  
 - Component has received forces resulting from a fall of factor 1 or more.  
 - Component has had more than 6 months of intensive use, 12 months of normal use, 10 years of occasional use.

The inspector accepts no responsibility in the case of omission or inaccuracy in the information concerning the checking of the component's history, which must be done by the client.

SCAN FOR INSTRUCTION

Model \_\_\_\_\_  
 Colour \_\_\_\_\_  
 Serial nr \_\_\_\_\_

Stretcher meets partially the requirements of DIN EN 1865:1995, DIN EN 1498:2007 and DIN 354:2002

[www.skedco.eu](http://www.skedco.eu)

Colour \_\_\_\_\_  
 Individual Serial Number \_\_\_\_\_

For Example the Serial Number is 1745012

17= year of manufacturing / 45=week off purchase / 012=unique number

Colour: OR=International Orange / CB=Coyote Brown / GR=Green

**By scanning the QR Code with a free app on your smartphone, you can instantly download the manual on your phone or tablet.**

TRANSLATION

(1) **Type Examination Certificate**

(2) No. of the Type Examination Certificate: **ZP/B061/16** replaces ZP/B029/11 PZ

(3) Product: **Rescue System**  
 Type: **Sked Rescue System**

(4) Manufacturer: **Felix Trading BV**

(5) Address: **Industrieweg 36-04, 3606 AS MAARSSSEN, NETHERLAND**

(6) The design of this product and any acceptable variation thereto are specified in the schedule to this type examination certificate.

(7) The certification body of DEKRA EXAM GmbH certifies that this product comply with the fundamental requirements of the standard listed under item: 8 below. The examination and test results in the test and assessment report PB 16-104.

(8) The examined types meet partially the according requirements of

**DIN EN 1865-1:2010      DIN EN 1498:2007      DIN EN 354:2010**

(9) This Type Examination Certificate relates only to the design, examination and tests of the specified product in accordance to the standard list. Further requirements apply to the manufacturing process and supply of this personal protective equipment. These are not covered by this certificate.

(10) This Type Test Certificate is valid until 2021-05-30.

DEKRA EXAM GmbH  
 Bochum, 2016-05-31

signed: Wiegand  
 Certification body

signed: Stickdom  
 Special services unit

We confirm the correctness of the translation from the German original.  
 In the case of arbitration only the German wording shall be valid and binding.

*Wiegand*  
 Certification body

*Stickdom*  
 Special services unit

---

# SKED<sup>®</sup>

## THE COMPLETE RESCUE SYSTEM

Felix Safety BV  
Industrieweg 36 - 04  
3606 AS Maarsse  
The Netherlands

T+ 31 346 - 555 622  
F+ 31 346 - 555 760

[info@felixsafety.com](mailto:info@felixsafety.com)  
[www.felixsafety.com](http://www.felixsafety.com)

---